


# HFIVD<sup>R</sup> Extraction & Detection Kit Operation flow chart

## 1.Extraction kit



- ① Put the sample into a clean lab blender bags.
- ② Add NS (1-3 times the volume of the sample) to the sample and grind sufficiently.
- ③ Pipet 400  $\mu$ L ~1.5mL homogenate in a clean 1.5 mL microcentrifuge tube, and centrifuge for 15~20s at 6,000~12,000rpm (when use constant speed centrifuge) or 1min at 2,000rpm (when use adjustable speed centrifuge).
- ④ Pipet 300  $\mu$ L supernatant in a clean 1.5 mL microcentrifuge tube.

- ① Add 450  $\mu$ L Lysis Buffer to the supernatant, mix by vortex mixer and incubate at room temperature for 5min.
- ② Centrifuge the mixture for 1 min at 12,000rpm


- ① To bind DNA/RNA, apply the mixture to the HFquick column, and centrifuge for 1min at 12,000rpm .
- ② Discard flow-through and place HFquick column back in the same collection tube.

- ① Add 700  $\mu$ L of Wash Buffer to HFquick column and centrifuge for 1min at 12,000rpm .
- ② Discard the flow-through and centrifuge the HFquick column for an additional 3min at 12,000rpm .

Place HFquick column into a clean 1.5 ml microcentrifuge tube.

- Add 30  $\mu$ L Buffer EB to the center of the HFquick membrane, let the column stand for 3 min, and then centrifuge for 1 min at 12,000rpm. Store DNA/RNA .
- Note:
- ✓ For DNA: stored at  $-20^{\circ}\text{C}$ .
  - ✓ For RNA: stored at  $-70^{\circ}\text{C}$  or  $-20^{\circ}\text{C}$  when less than 24h.

## 2.Detection kit



Preheat the thermocycler : set the PCR program .

Add 50 $\mu$ L deionized water to the Freeze-dried Positive control DNA/RNA tube .

Add 22.5 $\mu$ L Reagent B to the Reagent A tube .

- ① Add 2.5 $\mu$ L negative control, positive control, and sample separately.
- ② Cover the tube.

- ① Vortexing or finger flicking, then centrifuging for 30s-1min at  $\geq 12,000$ rpm to get rid of the bubble .
- ② Finally transfer the tubes to preheated thermocycler and start the program .